

Seating Accessories

- **Material:**
Memory foam
- **Included Hardware:**
3 Bolts
- **Dimensions:**
Small - 9" wide
Large - 11" wide

Soft Headrest Pad

Curved



DESCRIPTION	ITEM #
Small.....	053508
Large.....	053507

- **Description:**
Multi-axis headrest positioning hardware

- **Included Hardware:**
Mounting clamp

Headrest Hardware



DESCRIPTION	ITEM #
Headrest Hardware.....	053501
Replacement Clamp.....	053505

- **Description:**
Flip down, height adjustable

• **Dimensions:**

	width	depth	height
small	2 1/4"	4 1/4"	3"
medium	3 1/4"	5 1/4"	3"

- **Included Hardware:**
2 Bolts
2 Washers

Abductor



DESCRIPTION	ITEM #
Small.....	05354
Medium.....	05355

- **Material:**
Birch wood

- **Dimensions:**
8" Wide

- **Disclaimer:**
Transfer Boards are designed to help transfer to and from a wheelchair.
Weight capacity 250 lbs.

Transfer Board



DESCRIPTION	ITEM #
	Hand Holes No Hand Holes
Tapered, 26".....	847102005935
Tapered, 34".....	847102005386
Non-Tapered, 35".....	056102
Non-Tapered, 40".....	847102019130

Armrest Pad



DESCRIPTION	ITEM #
-------------	--------

Full Length.....	847102015279
Desk Length.....	847102015286

Seating Accessories

- **Description:**
Sold individually

Torso Pad



DESCRIPTION	ITEM #
-------------	--------

Small.....	05261002
Medium.....	05261102
Large.....	05261202

- **Description**
Sold individually

- **Dimensions:**
Small - 2 1/2" x 4 1/2"
Medium - 3" x 6"
Large - 4" x 8"

- **Included Hardware:**
2 Bolts

Shoe Holders

Standard



DESCRIPTION	ITEM #
-------------	--------

Small.....	0526104
Medium.....	05262
Large.....	05263
Replacement Straps	052644

- **Description:**
Sold as a pair

- **Dimensions:**

	length	heel width	toe width
small	5 1/4"	3"	3 1/4"
medium	7 1/4"	3 1/4"	3 3/4"
large	9 1/4"	4"	4 1/2"

- **Included Hardware:**
4 Bolts
4 Nuts

Oxygen Cylinder Holder



DESCRIPTION	ITEM #
-------------	--------

Oxygen Cylinder Holder.....	06091
-----------------------------	-------

- **Description:**
Fits most standard wheelchairs. Designed to hold D and E cylinders. Connects to the wheelchair's anti-tipper.

- **Disclaimer:**
Oxygen Cylinder Holder is not designed as a stabilizing device. It's primary function is to hold an oxygen cylinder.